Tissue, Section Disease, Hypertension, Kidney Tissue Section

PRODUCT OVERVIEW

Product Name	Tissue, Section Disease, Hypertension, Kidney, Tissue Section
Full Product Name	Tissue, Section, Human Disease, Hypertension, Kidney (Frozen)
Purity/Purification	Molecular Biology Grade
Preparation and Storage	Store at -70°C.
Donor Information	86 year-old, adult female kidney. Clinical Identification: Hypertension.
Product Categories/Family for Tissue, Section Disease, Hypertension, Kidney tissue section	Cells; Lysates; Plasma; Serum; Tissues; Tissues-Sections
Related Product Information for Tissue, Section Disease, Hypertension, Kidney tissue section	Tissues were snap frozen in liquid nitrogen immediately after excision, and embedded in OCT. All frozen tissue products, including frozen tissue arrays, frozen tissue sections, and frozen tissue panels, can be used for both immunohistochemistry and in situ hybridization. These products are ideal for rapid cellular localization of RNA and protein expression, especially for antibodies that don't recognize paraffin embedded tissue. At least one of the tissue slides from each lot was stained with H & E to ensure the quality. The tissue sections are 5-10um in thickness, mounted on positively charged glass slides, and fixed by cold acetone.
Other Notes	Small volumes of Tissue, Section Disease, Hypertension, Kidney tissue section vial(s) may occasionally become entrapped in the seal of the product vial during shipment and storage. If necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the container`s cap. Certain products may require to ship with dry ice and additional dry ice fee may apply.
Searchable Terms for Tissue, Section Disease, Hypertension, Kidney purchase	MBS640609 is a tissue section targeting Tissue, Section Disease, Hypertension, Kidney. Also known as Tissue, Section, Human Disease, Hypertension, Kidney (Frozen).



Testing Data of Tissue, Section Disease, Hypertension, Kidney tissue section
Human

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

MyBioSource | Tel: +1.858.633.0165 | Fax: +1.858.633.0166 | mybiosource.com/contact-us
Online datasheet: https://www.mybiosource.com/tissue-section/tissue-section-disease-hypertension-kidney/640609